

New concept proposal

Library Preparation

Author	Eelke van der Horst, Femke Kopmels	Date last updated	31/10/2023
Project	General Interest	Contact	DCC
Dataset release	2024.1	Consulted expert	-

1 Rationale



Library Preparation is a type of *Sample Processing* that is essential to prepare a DNA or RNA sample for sequencing using a next generation sequencing assay. It usually consists of isolating nucleic acids from the input sample, fragmenting the targeted sequences, attaching adapter sequences, and optional amplification, resulting in a sample of DNA fragments i.e. the NGS library. This library is then sequenced using the sequencing instrument of the sequencing assay. Participants that characterise this process are the Library preparation kit and Target enrichment kit, i.e. ready-to-use kits consisting of reagents and a fixed protocol to perform a step in library preparation. The choice of these kits depend on the intended type of NGS experiment, i.e. whole genome sequencing, whole exome sequencing, or targeted sequencing, and the sequencing platform that will perform the subsequent sequencing assay. Other essential attributes are the intended insert size (for paired-end sequencing) and whether PCR amplification was used.

2 Comparison to other standards/data models

2.1 FAIR Genomes

FAIR Genomes has a 'Sample preparation' module which describes the same experimental step as *Library Preparation*; a meaning binding between the 'fg:Sample_preparation' class and the *Library Preparation* concept is therefore proposed. While the label 'Sample preparation' in the context of FAIR Genomes can only mean preparation of a sample for nucleic acid sequencing, this label is ambiguous in the context of the SPHN, since it can also mean preparation for other (omics) assays. We therefore denote this concept with a more specific label 'Library Preparation'. The elements for the library preparation kit and target enrichment kit from FAIR Genomes are mirrored in this concept, reusing the value sets as defined by FAIR Genomes. The 'intended insert size' is also mirrored from FAIR Genomes.

Note: FAIR Genomes is referred to in the SPHN meaning bindings with the prefix 'fg'.

A project of	 <p>SAMWASSM Schweizerische Akademie der Medizinischen Wissenschaften Académie Suisse des Sciences Médicales Accademia Svizzera delle Scienze Mediche Swiss Academy of Medical Sciences</p>	 <p>SIB Swiss Institute of Bioinformatics</p>	<p>SIB Swiss Institute of Bioinformatics PHI Personalized Health Informatics Group www.sphn.ch dcc@sib.swiss</p>
--------------	--	--	--

2.2 ENA

The ENA ‘Experiment’ object has fields that describe the library preparation, including insert size, library name, library source, library selection, library strategy, and an optional free text library description. However, information on library source is part of the upstream sample processing and specimen, while the library strategy is reflected in the type of sequencing assay that this library preparation is part of. Any special library selection is either implicit to the NGS kits used or should be indicated with a preceding sample processing step.

2.3 OBI and EFO

The Ontology of Biomedical Investigations (OBI) has a ‘sample preparation for sequencing assay’ (OBI:0001902) concept, defined as “A sample preparation for assay that preparation of nucleic acids for a sequencing assay”. This corresponds to *Library Preparation* in case of NGS sequencing, however, the OBI class might also be used for non-NGS-based sequencing. It also has the ‘library preparation’ class (OBI:0000711). EFO reuses both OBI classes, however, defines more specific subclasses for library preparation than OBI, for instance ‘cDNA library construction’ and ‘single cell library construction’ (and more specific ‘CITE-seq’).

3 Concept information

Concept or concept compositions or inherited	General concept name	General description	Contextualized concept name	Contextualized description	Type	Standard	Value set or subset	Meaning binding	Cardinality for composedOf
concept	Library Preparation	process which results in the creation of a library from fragments of DNA	Library Preparation	process which results in the creation of a library from fragments of DNA	Sample Processing			OBI:0000711 library preparation	
composed Of	kit code	coded information specifying the kit associated to the concept	library preparation kit code	pre-filled, ready-to-use reagent cartridges intended to improve chemistry, cluster density and read length as well as improve quality (Q) scores for this sample. Reagent components are encoded to interact with the sequencing system to validate compatibility with user-defined applications.	Code	EFO, GENEPIO, FAIRGenomes or other		fg:Sample_preparation_Library_preparation_kit	0:1
composed Of	target enrichment kit code	coded information specifying the target enrichment kit	target enrichment kit code	indicates which target enrichment kit was used to prepare this sample. Target enrichment is a pre-sequencing DNA preparation step where DNA sequences are either directly amplified (amplicon or multiplex PCR-based) or captured (hybrid capture-based) in	Code	EFO, GENEPIO, FAIRGenomes or other		fg:Sample_preparation_Target_enrichment_kit	0:1

		associated to the concept		order to only focus on specific regions of a genome or DNA sample.					
composed Of	intended insert size	intended insert size associated to the concept	intended insert size	in paired-end sequencing, the DNA between the adapter sequences is the insert. The length of this sequence is known as the insert size, not to be confused with the inner distance between reads. So, fragment length equals read adapter length (2x) plus insert size, and insert size equals read length (2x) plus inner distance.	Quantity				0:1
composed Of	gene panel	gene panel associated to the concept	gene panel	collection of genes that are the focus of sequencing	Gene Panel				0:1
inherited	code	coded information specifying the concept	code	code specifying the type of library preparation	Code	OBI, EFO or other	for OBI: descendant of: OBI:0000711 library preparation		0:1
inherited	standard operating procedure	standard operating procedure associated to the concept	standard operating procedure	standard operating procedure that was followed for this library preparation	Standard Operating Procedure				0:1
inherited	input	input associated to the concept	input	the sample for which a library is created	Sample				0:n

inherited	output	output associated to the concept	output	the NGS library that is produced	Sample				0:1
inherited	start datetime	datetime at which the concept started	start datetime	start of library preparation	temporal				0:1
inherited	predecessor	a preceding process associated to the concept	predecessor	process preceding this library preparation	Sample Processing			RO:0002087 immediately preceded by	0:n
inherited	quality control metric	quality control metric associated to the concept	quality control metric	quality control metric related to the output of the library preparation	Quality Control Metric			EDAM:data_3914 quality control report	0:n

General concept name	Cardinality for concept to Administrative Case	Cardinality for concept to Data Provider	Cardinality for concept to Subject Pseudo Identifier	Cardinality for concept to Source System
Library Preparation	0:n	1:1	0:n	1:n

4 Impact on the SPHN Dataset

-

5 Discussion

Library preparation can be considered an integral part of performing an NGS sequencing experiment. However, since the configuration/parameters of these experimental processes can vary separately, and also because they can be executed at different facilities at different times by different people, while also having their independent quality control metrics, we represent these as separate concepts that are part of the *Sequencing Assay*. Note that the parent *Sequencing Assay* and its *Sequencing Instrument* are tightly bound to this concept, because these influence the possible choices for the kits used for library preparation and constrains the intended read length and insert size.

An alternative naming that was considered is 'sample preparation for sequencing assay'. However, in the context of Omics/NGS sequencing, library preparation seemed more specific.

Note that as with other (experimental) processes, individual sub-steps can be provided, for instance for DNA extraction or amplification.

Since many of the composedOf concepts might not be known, for instance which kit was used, these have minimum cardinality 0.

6 Example

Library preparation as part of a WGS sequencing assay

preparation kit code: **Illumina TruSeq DNA PCR-Free**

start datetime: **2023-06-28**

intended insert size

value: **350**

unit: **{#}**

input:

identifier: **DNASample1**

collection datetime: **2023-06-26**

material type code: **258566005 [Deoxyribonucleic acid specimen]**

output:

identifier: **Lib1**

collection datetime: **2023-06-28**

material type code: **258566005 [Deoxyribonucleic acid specimen]**

predecessor:

code: **104166004 [Nucleic acid molecular isolation or extraction method (procedure)]**

input:

A project of	 <p>Schweizerische Akademie der Medizinischen Wissenschaften Académie Suisse des Sciences Médicales Accademia Svizzera delle Scienze Mediche Swiss Academy of Medical Sciences</p>	 <p>Swiss Institute of Bioinformatics</p>	<p>SIB Swiss Institute of Bioinformatics PHI Personalized Health Informatics Group www.sphn.ch dcc@sib.swiss</p>
--------------	---	--	--

identifier: **BloodSample1**
 collection datetime: **2023-06-26**
 material type code: **19297000 |Blood specimen|**
 output:
 identifier: **DNASample1**
 collection datetime: **2023-06-26**
 material type code: **258566005 |Deoxyribonucleic acid specimen (specimen)|**
 standard operating procedure:
 name: **DNA isolation from blood sample**
 description: **To describe how DNA from blood must be isolated and stored.**
 version: **V.1.1**
 data file:
 format code: **EDAM: format_3508 |PDF|**
 code: **EXAMPLE01**